Payvand

Clinical Specialty Lab.

Hematology Workshop

Quality Assurance with emphasis on IQC of Cell Counters Day 2

Behzad Poopak, DCLS PhD

Head Payvand Clinical & Specialty Lab.

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Purpose & Method of IQC

- Once a performance baseline has been established, its control is based on methods that will **detect loss of accuracy** or **loss of precision**, or both.
- Control methods must be sensitive enough to reveal loss of performance that could compromise patient assay values, but they ***cannot be oversensitive*** to the point of signaling error when none exists.

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- **As a minimum**, a performance control system must provide the following information:
 1. **Analytical imprecision** has not deteriorated with time or use.
 2. **Bias** has remained within the limits established at the time of calibration.
- The term "**internal**" indicates quality control procedures that take place within the laboratory immediately **before**, **during**, or immediately **after** an analysis run.

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- The procedures must be carried out as an *integral part of patient assays*.
- The definition of a "*run*" should be related to **work loading** and to the **stability (drift)** of the analyzer, not simply to elapsed time.
- Control assays performed **only once daily** may not provide enough information to ensure accuracy and consistency of results.
- *When the analyzer is restarted after complete shut-down*, **precision** and **set point** should be verified to ensure that satisfactory analytical capability has not been lost.

ERRORS: Random and Systematic

- *All measurements performed in a hematology laboratory contain inherent variability*
- Random errors – introduce *variance*
 - » increase the scatter of values about the true value
 - » result of chance (eg. sampling errors)
 - » do not affect an entire batch of specimens
 - » are not be detected by control samples
- Systematic errors - introduce *bias*
 - » not due to chance
 - » eg. deteriorating reagents or loss of calibration

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In 1931,

Dr. Walter Shewhart, a scientist at the Bell Telephone Laboratories, proposed applying **statistical based control charts** to interpret industrial manufacturing processes.

Dr. Walter A. Shewhart

In 1950,

**S. Levey &
E.R. Jennings**
suggested the use in
the clinical laboratory.

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
The chart will suggest that there is a fault in technique, instrument, pipette or reagents if one of the following occurs :

First check that the material itself has not become infected or in other ways has begun to deteriorate

Levey-Jennings Control chart

In control: values oscillate around mean, within ± 2 SD

Test Replication

- The mean of replicate assays of control material is more informative than a single measurement.
- **Paired assays** are recommended. Not only does the mean of the pair give a **more exact value**, but the difference between members of a pair also gives an **estimate of imprecision**.
- **Overlap** is a special case of test replication. When a **new lot of control material is received**, it should be tested in parallel with the currently used lot **for four days**.
- This will establish the ratio of the assigned values of the incoming lot to the assigned values of the currently used material.
- **Differences in assigned values of the two lots should not be taken as an indication to recalibrate.**

When a control point is outlier...

- You can use two types of rules for interpretation:

WHO & Westgard

Warning rule = use other rules to inspect the control points

Rejection rule = “out of control”

- Stop testing
- Identify and correct problem
- Repeat testing on patient samples and controls
- **Do not report patient results until problem is solved and controls indicate proper performance**

Westgard Rules

- Multi-rule QC
- It uses a combination of decision criteria, usually 5 different control rules to judge the acceptability of an analytical run.

Explanation of Individual Rules

- **1_{2s}** : One control measurement exceeding 2 standard deviations of control limits either above or below the mean. This rule is used as a **warning rule** to trigger careful inspection of the control data.
- **1_{3s}** : This rule is commonly used with a Levey-Jennings chart when the control limits are set as the mean ± 3 standard deviations of control limits. A run is **rejected** when a single control measurement exceeds the mean ± 3 control limits.
- **2_{2s}** : The control run is **rejected** with 2 consecutive control measurements 2 standard deviations of control limits on the same side of mean with this rule.
- **R_{4s}** : This rule **rejects** a run if two control measurements in a group exceed the mean with a 4 standard deviation difference between the 2 controls
- **4_{1s}** : This rule **rejects** a run with the 4th consecutive control measurement exceeding 1 standard deviation on the same side of the mean.
- **10_x** : This **rule rejects** a control run when there are 10 consecutive controls on the same side of the mean.


WHO Rules

1 control value outside the mean +/- 2SD	Warning
1 control value outside mean +/- 3SD	Reject : SE or RE
2 consecutive controls exceed mean +/- 2SD	Reject : SE
4 consecutive controls exceed mean + SD or mean - SD	Reject : SE
6 consecutive controls on one side of the mean	Warning : SE

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
In control: values oscillate around mean, within +/- 2 SD


Out of control : values drift, problem progressively developing

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
Out of control, shift: abrupt change, values oscillate around new mean

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Out of control : Increased dispersion, random errors eg. technique

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Moving Averages (*Bull's method*), Patients Data

- large laboratories (at least 100 blood count / day)
- Assumes the population sampled each day remains constant
- Therefore the calculated indices (MCV, MCH and MCHC) remain stable
- As **any significant change** may indicate : *instrument or technical fault*
- As accurate as 4C preparations
- For manual method use MCHC

Moving Averages

- Before setting up this procedure it is necessary to establish the SD.
- This is done by calculating the **mean on 11** consecutive working days & calculating the SD

$$SD = \sqrt{\frac{\sum(x - \bar{x})^2}{n - 1}} \text{ where } x = \text{daily mean, } \bar{x} = \text{mean of daily means and}$$


$n =$ number of days. In this case, therefore, $n - 1 = 10$.

- Plot the data on a graph
- The mean rbc indices is calculated at the end of each day
- Determine mean indices values for each batch of 20 patients, plot on control chart


Moving Averages

- Results are valid only if the population does not significantly from *day to day*
- Tests are **not selectively biased** (specific outpatient clinics) (exclusion or median)
- Any batch of 20, no more than 7 should come from one clinical source or have the same clinical condition.
- If the test is being performed satisfactorily, the mean will not vary by more than 2SD on any day.
- The acceptable range is **+/- 3%** of the mean.
- ***Moving averages cannot be used to monitor leukocyte or platelet counts, as those counts are variable within a patient population.***

Moving Averages

Moving averages - low Hemoglobin

Changes in the moving averages graphs indicate where the problem might be in the system. eg. If the light source for Hb is becoming weak, then the calculated MCH and MCHC values will fall

Moving averages

<i>CAUSE</i>	MCV	MCH	MCHC
Low Hb	no change	low	low
High Hb	no change	high	high
Low RCC	high	high	no change
High RCC	low	low	no change
Low Hct	low	no change	high
High Hct	high	no change	low